

Structomer Hard Splint

Instructions for Use (IFU)

3D Printing Resin for Hard Splint

For Professional Dental Use Only

1. Product Description

Structomer Hard Splint is a Class II, biocompatible 3D printing resin formulated for the fabrication of full maxillary arch custom bite splints and corrective oral appliances. It is designed for use with CAD/CAM systems and validated 3D printing and post-curing devices.

2. Composition

Structomer Hard Splint resin consists of:

- Dimethacrylate-based resin
 - Photoinitiator
 - Monomer (Diluents)
 - Additives
-

3. Intended Users

This product is intended for use exclusively by trained dental professionals, including licensed dentists and dental laboratory technicians.

Sales are restricted to:

- Licensed dental professionals
- Dental supply distributors
- Educational institutions
- Government dental facilities

Note: This product is labelled for use by licensed practitioners only, unless permitted by local or state law.

4. Indications for Use

Structomer Hard Splint resin is indicated for:

- Custom full-arch bite splints
- Mouthguards for bruxism (teeth grinding)
- Anti-snoring oral appliances

This resin is to be used in conjunction with a complete CAD/CAM system, which includes:

- Oral impression or intraoral scan
- Digital design file
- 3D additive manufacturing printer (stereolithography)
- Post-curing equipment

Use only with the validated printers and post-curing devices listed in this guide, following all manufacturer recommendations.

5. Contraindications

- Do not design bite splints with thicknesses below **0.75 mm**, as this may result in breakage or deformation during removal from the patient's mouth.

6. Design Orientation & Support Recommendations

6.1 Orient bite splints with the adaptive (occlusal) surface angled **downward at 45°**.

6.2 Place support structures around the **perimeter of the ridge**. Avoid placing supports on adaptive or gingival contact areas.

6.3 **Minimum support diameter:** 0.3 mm at contact point

6.4 **Minimum support height:** 2 mm (to ensure easy removal from build plate)

7. Mixing Instructions

Thorough mixing is essential due to the different weights of resin components.

- **7.1 Resin in Vat:** Stir gently with a silicone spatula if idle for more than 4 hours.
- **7.2 From Bottle:** Shake or stir with a plastic spatula for at least 1 minute, or use a bottle roller for 10 minutes.
- **7.3 Long-term Storage:** If the resin has been stored for over 1 month, roll the bottle for **1 hour** before use.

8. Post-Processing Instructions

8.1 Remove build plate after printing.

8.2 Clean parts in **98% IPA**:

- 1st bath (ultrasonic): 2 minutes
- 2nd bath with fresh IPA: 2 minutes
Do not exceed 4 minutes total IPA exposure to avoid weakening the material.

8.3 Use compressed air to dry parts and remove uncured resin.

8.4 Repeat cleaning until a **matte, shine-free** finish is achieved.

8.5 Post-cure using a validated device according to recommended time and temperature (see Section 10).

Post-curing is required for FDA compliance.

9. Validated 3D Printers

Use only the following validated printers (wavelength 405nm):

- Structo Dentaform

10. Validated Post-Curing Devices

- **Formlabs Form Cure, Fisrt Gen.:** 5mins at 60 °C
 - Light Source: 13 multi-directional LEDs
 - LED Power: 39 W
 - LED Radiant: 9.1 W
 - LED Wavelength: 405 nm

- Structo Postcuring Units: 15mins

Light source: Six 18W/71 lamps (Dulux L Blue) and six 18W/78 lamps (Dulux blue UV-A)

No vacuum or inert gas environment is required.

11. Chairside Adjustments & Impressions

11.1 Adjustments:

Use standard dental carbide burs to modify the bite splint. Smooth sharp areas using polishing wheels.

11.2 Impression Taking:

Follow the impression material manufacturer's instructions and ensure full set before removal.

12. Storage and Handling Conditions

12.1 Light Sensitivity:

Keep resin away from direct sunlight and ambient light. Cover resin vats when not in use.

12.2 Temperature:

- **Storage:** 65°F – 85°F (18°C – 29°C)
 - **Printing:** Ideally at 30°C
If stored in cold environments, warm the sealed resin bottle in a water bath before use.
-

13. Long-Term Storage & End-of-Day Procedures

- Return unused resin to the original bottle using a **fine mesh filter**.
 - Store in **original packaging** at room temperature in a **cool, dry, dark place**.
 - Keep container tightly closed when not in use.
 - **Shelf Life:**
 - In bottle: Up to **12 months**
 - In printer: Up to **1 week** with hood closed
 - Do **not expose to UV light**.
 - Standard transport conditions apply; **no special hazardous material protocols** required.
-

14. Disposal Guidelines

Dispose of resin in accordance with **federal, state, and local** regulations.

- Refer to Structo soft splint MSDS for hazardous waste classification.
 - **Cure liquid resin** before disposal: pour into a clear container and expose to sunlight until hardened.
 - Cured resin may be discarded in regular trash.
-

15. Legal Notice

Structo assumes no liability for misuse of this product, including:

- Use of non-validated equipment
 - Deviation from the instructions provided
Such actions may affect the performance, safety, or effectiveness of the final restoration.
-

16. Packaging Information

- **Available Variant:**
 - **Color:** Clear
 - **Unit Size:** 1kg & 5kg bottle
-